

VISTAZO™

Hydrophilic Lens
Quadra Haptic
Aspheric
Clear / Yellow



Highly biocompatible
Square Edge Design to minimize PCO
Excellent optical performance
Excellent stability and centration due to the unique design of the haptic

OPHTHALMIC
Solution, Inc.
Manufacturer of Intraocular Lenses

- The special aspheric design ensures exceptional optical performance, even under tilt or decentration.
- Optimized to accommodate various corneal aberration profiles.

Technical Specifications

Model	QL360110CR6A	QL360110EY6A
Overall Diameter	11.00 mm	11.00 mm
Optic Diameter	6.00 mm	6.00 mm
Haptic Angle	5°	5°
Haptic Design	Quadra Haptic	Quadra Haptic
Optic Type	Biconvex, Aspheric	Biconvex, Aspheric
Material	Clear Hydrophilic Acrylic	Extended Yellow Hydrophilic Acrylic
UV Blocker	Yes	Yes
Blue light filter	No	Yes
Refractive Index	1.46	1.46
Diopter Range	Monofocal >5.00 upto 45.00 (0.5 increment)	Monofocal >5.00 upto 45.00 (0.5 increment)

Optimized IOL constants for the IOL Master QL360110CR6A / QL360110EY6A

Nominal	118.4
Haigis	a0 = 1.36 a1 = 0.40 a2 = 0.10
Hoffer Q	pACD = 5.54
Holl. 1	sf = 1.74
SRK / T	A = 118.8
SRKII	A = 119.2

